

Clonidine Cleanup

Use $\frac{3}{4}$ - 1 tablet and crush into a powdered form.

Place the powdered sample into a small labeled test tube and dissolve with $\frac{1}{2}$ - 1mL of hot water.

Vortex and let stand for $\frac{1}{2}$ hr.

Remove and filter the water and place into another small labeled test tube and then add sodium bicarbonate.

Mix the solution and perform a micro-extraction with $\frac{1}{2}$ - 1mL of chloroform.

Place the solvent into a residue vial for GC/MS confirmation.

Annie